

ENVIRONMENTAL ANALYSIS ASSOCIATES, INC.

EXAMPLE REPORT EXCERPTS

Automated Scanning Electron Microscopy / SDD X-ray dust analysis



The excerpts on the following pages are from 3 different types of dust analysis reports generated by the new Automated SEM procedure developed by Environmental Analysis. This analysis capability is now located in both the Bay City Michigan, and San Diego California laboratory locations. The particle size/mass distribution, and chemistry of complex airborne and surface dust samples can be simultaneously quantified, and then "classified" into a profile of probable emission sources at an economical price (\$250 - \$350/sample).

Example Applications :

- Identification of IAQ dust contamination sources
- Wildfire and structure fire particle analysis
- Respirable quartz in airborne coal dust
- Inhalable / respirable particle size analysis
- Monitoring of hospital OR's and clean rooms
- Product / dust comparison analysis

Particle Classification	Particles Counted	Mean size (um)	* Specific Gravity	Numerical %	* Mass Ratio %
Ca oxide	61	13.7	3.00	30.5%	60.9%
Ca carbonate	49	8.3	3.00	24.5%	16.0%
CaMgSiK carbonate	44	4.1	3.00	22.0%	1.2%
MgK carbonate	3	10.2	3.00	1.5%	1.3%
CaMgK phosphate	13	6.6	3.00	6.5%	2.6%
Si carbon	3	5.0	1.50	1.5%	0.1%
Quartz	3	5.7	2.60	1.5%	0.1%
Ca sulfate	3	3.8	3.00	1.5%	0.1%
Mixed Ca silicate	1	16.2	3.00	0.5%	1.0%
MgAlKCaFe silicate	2	20.4	3.00	1.0%	12.8%
Unclassified	18	7.1	2.00	9.0%	4.1%

Wildfire fire ash – Particle classification summary table



For more information contact Daniel Baxter at Environmental Analysis Associates, Inc.

San Diego, California

5290 Soledad Road
San Diego, CA 92109

Bay City, Michigan

306 5th Street, Suite 400
Bay City, MI 48708

Phone: 858-272-7747

Email: dbaxter@eaalab.com

Website: eaalab.com

Automated Scanning Electron Microscopy - Dust Analysis Photo Report

Client Name : Research

Contact : Daniel Baxter

Client Project# : Fire pit ash 4-14-15

Client Sample # : 041015-1


Sample Description : Oak fire ash sample from Michigan fire pit

Analysis Method : Automated SEM/EDS


Analysis Date : 8/14/15

EAA Project # : Fire pit


EAA Sample # : Fire ash-1



204x - Secondary electron image of analyzed field



204x - Backscatter image of analyzed field



Calcium "oxide" crystals - 490x



Calcium "oxide" crystals - 2010x

Automated Scanning Electron Microscopy - Dust Analysis Report

Summary Page

Client Name : Research	Analysis Date : 8/14/15
Contact : Daniel Baxter	EAA Project # : Fire pit
Client Project# : Fire pit ash 4-14-15	EAA Sample # : Fire ash-1
Client Sample # : 041015-1	
Sample Description : Oak fire ash sample from Michigan fire pit	Fields Counted : 2
Analysis Method : Automated SEM/EDS	Field area cted (mm ²) : 0.095
Analysis Magnification : 204	Particles / mm ² : 2103
Min. size threshold (um) : 2.0	Scale (um/div.) : 1
Max. size threshold (um) : 100	Total particles counted : 200

SUMMARY CONCLUSIONS

The fire ash was fully "decarbonized" at high temperature and is primarily composed of Calcium oxide / hydroxide, and Calcium carbonate and Magnesium, Potassium (K) salts. Minor amounts of Calcium/Magnesium silicates and phosphorus (P) oxides are also present. The Calcium oxide / oxalate crystals are the most common indicator particles for wood "ash".

Particle Classification	Particles Counted	Mean size (um)	* Specific Gravity	Numerical %	* Mass Ratio %
Carbon H	1	8.3	1.50	0.5%	0.1%
M carbon	9	4.0	1.50	4.5%	0.2%
Ca oxide	61	13.7	2.00	30.5%	59.7%
Ca carbonate	53	8.2	2.00	26.5%	16.6%
CaMgSiK carbonate	46	4.2	2.00	23.0%	1.4%
MgK carbonate	3	10.2	2.00	1.5%	1.2%
CaMgK phosphate	13	6.6	2.00	6.5%	2.5%
Quartz	4	7.9	2.00	2.0%	0.8%
Ca sulfate	3	3.8	2.00	1.5%	0.1%
Mixed Ca silicate	1	16.2	2.00	0.5%	1.0%
MgAlKCaFe silicate	2	20.4	2.00	1.0%	12.5%
Unclassified	4	10.8	2.00	2.0%	3.8%

* Mass ratios are based on the average particle size & theoretical specific gravity. If the specific gravity is unknown and not "carbonaceous", it will automatically be assigned a density value of 2.0
 "Unclassified" particles fall outside of the established chi square rule fit criteria


Analysis Method : Automated EDAX Genesis X-ray particle analysis software integrated with a data libraries developed by Environmental Analysis

Analyst : Daniel M. Baxter

Date : 8/14/15

Analysis Date : 8/14/15
EAA Project # : Fire pit
EAA Sample # : Fire ash-1

Chemical Component	Percentage
Ca oxide	59.7%
Ca carbonate	16.6%
MgAlKCaFe silicate	12.5%
Unclassified	3.8%
Carbon H	0.1%
M carbon	0.2%
Ca sulfate	0.1%
Quartz	0.8%
CaMgK phosphate	2.5%
MgK carbonate	1.2%
CaMgSiK carbonate	1.4%
Mixed Ca silicate	1.0%



Automated Scanning Electron Microscopy - Dust Analysis Photo Report

Client Name : ABC Environmental

Contact : Mr. John Doe

Client Project# : 14-1000

Client Sample # : 110314-1


Sample Description : "Trace" white powder in HVAC system duct

Analysis Method : Automated SEM/EDS

Analysis Date : 8/14/15

EAA Project # : 14-0495

EAA Sample # : 1000-1



Field 1 107x - Flake-like corrosion particles

Automated Scanning Electron Microscopy - Dust Analysis Report

Summary Page

Client Name : ABC Environmental	Analysis Date : 8/14/15
Contact : Mr. John Doe	EAA Project # : 14-0495
Client Project# : 14-1000	EAA Sample # : 1000-1
Client Sample # : 110314-1	
Sample Description : "Trace" white powder in HVAC system duct	Fields Counted : 1
Analysis Method : Automated SEM/EDS	Field area cted (mm ²) : 0.127
Analysis Magnification : 107	Particles / mm ² : 315
Min. size threshold (um) : 3.0	Scale (µm/div.) : 1
Max. size threshold (um) : 1000	Total particles counted : 40

SUMMARY CONCLUSIONS

Sample is composed of a mixture of aluminum oxide and mixed aluminum, iron, and zinc corrosion oxide particles. These particles also contain low concentrations of chlorine / chlorides.

The dust is consistent with corrosion of the ducting and "galvanized" duct coatings.

Carbon H = Particle with greater than 60% carbon

Particle Classification	Particles Counted	Mean size (um)	* Specific Gravity	Numerical %	* Mass Ratio %
Carbon H	3	58.4	1.00	7.5%	2.9%
Al oxide	27	49.0	1.00	67.5%	61.1%
Zn oxide	2	50.9	1.00	5.0%	2.1%
M AlFeZn oxide	8	42.5	2.00	20.0%	33.9%
Unclassified					

* Mass ratios are based on the average particle size & theoretical specific gravity. If the specific gravity is unknown and not "carbonaceous", it will automatically be assigned a density value of 2.0
 "Unclassified" particles fall outside of the established chi square rule fit criteria

Analysis Method : Automated EDAX Genesis X-ray particle analysis software integrated with a data libraries developed by Environmental Analysis

Analyst : Daniel M. Baxter



Date : 8/14/15

Automated Scanning Electron Microscopy -Grapical Report - Mass & Size Distribution

Client Name : ABC Environmental
Contact : Mr. John Doe
Client Project# : 14-1000
Client Sample # : 110314-1
Sample Description : "Trace" white powder in HVAC system duct
Analysis Method : Automated SEM/EDS

Analysis Date : 8/14/15
EAA Project # : 14-0495
EAA Sample # : 1000-1

Estimated Mass %



Automated Scanning Electron Microscopy - Dust Analysis Photo Report

Client Name : Research

Contact : Dan Baxter

Client Project# : Mining rock

Client Sample # : Coal mine-1

Sample Description : Air sample - During operations

Analysis Method : Simulated 25mm PC filter analysis

Analysis Date : 7/22/15


EAA Project # : Research

EAA Sample # : 1000-1

Coal mine-1
Field 1
Magnification (X) 497



Coal mine-1
Field 2
Magnification (X) 497



Scanning Electron Microscopy Air Sample - Summary Report

Filter / Impaction Air Sample

Client Name : Research	Analysis Date : 7/22/15
Contact : Dan Baxter	EAA Project # : Research
Client Project# : Mining rock	EAA Sample # : 1000-1
Client Sample # : Coal mine-1	
Sample Description : Air sample - During operations	Fields / passes counted : 2
Filter media / size / type : Simulated 25mm PC filter analysis	Field area counted (mm ²) : 0.069
Analysis Magnification : 497	Collection area (mm ²) : 398.2
Scale (µm/div.) : 1	% of sample counted : 0.02%
Total particles counted : 342	Particles / mm ² : 4965
Sample volume (m ³) : 2.400	Particles / sample : 1976892
Min./Max. size range (um) : 3 / 165	Estimated Particles / m ³ : 823705

SUMMARY CONCLUSIONS - Major Constituents

Simulated air sample output using the bulk analysis of the Rom rock.
 Sample is primarily highly carbonaceous (Carbon H >30%), mixed Alumino-silicate minerals (Al silicate >50%), Quartz (~2%). Minor amounts of Calcium carbonate (Ca carbonate) and Magnesium Calcium Carbonate (MgCa carbonate) present.

Minor Constituents:

Particle Classification	# Cted	Mean (um)	Num. %	*Calc Mass %	*Spec Grav	Part. / Sample	Part / m ³	*Calc Mass ug/m ³	** Mean Mass ug/m ³
Carbon H	91	8.5	26.6%	22.6%	1.50	526015	219173	444.7	105.6
AlSi carbon	29	11.1	8.5%	19.1%	1.50	167631	69846	376.8	74.1
Ca carbonate	2	11.2	0.6%	0.9%	3.00	11561	4817	17.5	10.6
MgCa carbonate	1	5.9	0.3%	0.0%	3.00	5780	2408	0.8	0.8
Al silicate	106	8.8	31.0%	29.6%	2.00	612721	255300	583.1	183.0
AlK silicate	34	8.8	9.9%	6.6%	2.00	196533	81889	129.2	58.8
M silicate	42	10.4	12.3%	15.9%	2.00	242776	101157	313.9	118.6
Quartz	29	5.8	8.5%	1.9%	2.60	167631	69846	37.8	19.0
Unclassified	8	13.6	2.3%	3.4%	2.00	46243	19268	66.8	50.8
TOTALS	342					1976892	823705	1971	621

* The calculated mass/m³ is based on the sum total of each particle volume & theoretical specific gravity.

The ** Mean Mass is calculated using the Mean particle size for each category. Both calculations should be used as the likely range of estimated particle mass concentrations.

Analysis Method : Automated EDAX Genesis X-ray particle analysis

Analyst : Daniel M. Baxter

Date : 7/22/15

Analysis Date : 7/22/15
EAA Project # : Research
EAA Sample # : 1000-1

Mineral	Percentage
Al silicate	29.6%
Carbon H	22.6%
AlSi carbon	19.1%
M silicate	15.9%
ALK silicate	6.6%
Quartz	1.9%
Unclassified	3.4%
Ca carbonate	0.9%
MgCa carbonate	0.0%

